


Custom measurement form:
PCO ORTHOSES CARD

PATIENT ID _____

PRODUCT CODE _____

NOTES _____




! WARNING!

- Chest and upper limb measurements should be taken with patient lying down.
- Hip and lower limb measurements are taken with patient standing.
- All measurements should be taken in centimeters!
- If the finger put between the child's body and the PCO product meets resistance, the size of the product has been selected correctly.

IMPORTANT!

Only a correctly filled in measurement card is a guarantee of proper compression of the product. We are not responsible for wrongly taken measurements!

HOW TO MEASURE?



- I** Arms circumference at the level of the head of the humerus
- II** Shoulder circumference measured through the armpit and acromion
- III** Length measured from armpit to wrist
- IV** Height measured from the armpit to the greater trochanter of the femur
- V** Measurement from the manubrium of the sternum through the perineum to the cervical vertebra C7
- VI** Height measured from waist to ankle
- VII** Height measured from parineum to ankle
- VIII** Patient is height

REQUIRED MEASUREMENTS

PRODUCT	MEASUREMENT	LENGTH MEASUREMENT	CONTROL MEASUREMENT	ZIP
HEAD				
PCO-H-01	1			
TRUNK				
PCO-T-01	2, 3, 4, 5, 6, 14, 15, 16, 17, 18	F, G, H, I, J	I, II, IV, V, VIII	
PCO-T-02	2, 3, 4, 5, 6, 14, 15, 16, 17, 18	F, G, H, I, J	I, II, IV, V, VIII	
PCO-T-03	2, 3, 4, 5, 6, 14, 15, 16, 17, 18	F, G, H, I, J	I, II, IV, V, VIII	
PCO-T-04	2, 3, 4, 5, 6, 14, 15, 16, 17, 18	F, G, H, I, J	I, II, IV, V, VIII	
PCO-T-09	2, 3, 4, 5, 6, 14, 15, 16, 17, 18	F, G, H, I, J	I, II, IV, VIII	Y
PCO-T-10	2, 3, 4, 5, 6, 7, 8, 14, 15, 16, 17, 18	E, F, G, H, I, J	I, II, III, IV, VIII	Y
PCO-T-11	2, 3, 4, 5, 6, 7, 8, 9, 10, 14, 15, 16, 17, 18	C, D, E, F, G, H, I, J	I, II, III, IV, VIII	Y
PCO-T-16	2, 3, 4, 5, 6, 14, 15, 16, 17, 18	F, G, H, I, J	I, II, IV, VIII	Y
PCO-T-17	2, 3, 4, 5, 6, 7, 8, 14, 15, 16, 17, 18	E, F, G, H, I, J	I, II, III, IV, VIII	Y
PCO-T-18	2, 3, 4, 5, 6, 7, 8, 9, 10, 14, 15, 16, 17, 18	C, D, E, F, G, H, I, J	I, II, III, IV, VIII	Y
PCO-T-23	2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19	F, G, H, I, J	I, II, IV, V, VIII	Y
PCO-T-24	2, 3, 4, 5, 6, 7, 8, 14, 15, 16, 17, 18, 19	E, F, G, H, I, J	I, II, III, IV, V, VIII	Y
PCO-T-25	2, 3, 4, 5, 6, 7, 8, 9, 10, 14, 15, 16, 17, 18, 19	C, D, E, F, G, H, I, J	I, II, III, IV, V, VIII	Y
UPPER LIMB				
PCO-A-01	6, 7, 8, 9, 10, 11, 12, 13	A, B, C, D, E	III	
PCO-A-02	8, 9, 10, 11, 12, 13	A, B, C, D	III	
PCO-A-03	9, 10, 11, 12, 13	A, B, C	III	
PCO-A-07	6, 7, 8, 9, 10	C, D, E	III	
PCO-A-08	8, 9, 10	C, D	III	
LOWER LIMB				
PCO-L-01	16, 17, 18, 19, 20, 21, 22, 23, 24, 25	I, J, K, L, M	VI, VII, VIII	
PCO-L-02	16, 17, 18, 19, 20, 21, 22, 23, 24	I, J, K, L	VI, VII, VIII	
PCO-L-03	16, 17, 18, 19, 20, 21	I, J, K	VI, VII, VIII	
PCO-L-08	19, 20, 21, 22, 23, 24, 25	K, L, M	VII	
PCO-L-09	19, 20, 21, 22, 23, 24	K, L	VII	
PCO-L-10	19, 20, 21	K	VII	
PCO-L-11	23, 24, 25	L, M	VII	
PCO-L-12	24, 25, 26, 27	M, N, O, P	VII	

No.	MEASUREMENT	CIRCUMFERENCE	WIDTH	LENGTH
1	HEAD CIRCUMFERENCE			A
2	NECK CIRCUMFERENCE			B
3	NECK WIDTH			C
4	SHOULDERS WIDTH			D
5	SHOULDER-NECK WIDTH			E
6	ARM CIRCUMFERENCE			F
7	CIRCUMFERENCE BETWEEN 6 AND 8			G
8	ELBOW CIRCUMFERENCE WITH A SLIGHTLY BENT HAND			H
9	CIRCUMFERENCE BETWEEN 8 AND 10			I
10	WRIST CIRCUMFERENCE			J
11	THUMB CIRCUMFERENCE			K
12	HAND CIRCUMFERENCE			L
13	MCP JOINTS CIRCUMFERENCE			M
14	MEASUREMENT UNDER THE ARMPITS			N
15	MEASUREMENT BETWEEN 12 AND 14			O
16	WAIST CIRCUMFERENCE			P
17	HIP CIRCUMFERENCE			
18	MEASUREMENT AT THE PERINEUM HEIGHT			
19	GROIN CIRCUMFERENCE			
20	CIRCUMFERENCE BETWEEN 19 AND 20			
21	CIRCUMFERENCE 5 CM ABOVE THE KNEECAP			
22	KNEE CIRCUMFERENCE			
23	CIRCUMFERENCE 5 CM UNDER THE KNEECAP			
24	CALF CIRCUMFERENCE BETWEEN 22 AND 25			
25	ANKLE CIRCUMFERENCE			
26	HEEL CIRCUMFERENCE			
27	METATARSUS CIRCUMFERENCE			

CONTROL MEASUREMENT	
I	
II	
III	
IV	
V	
VI	
VII	
VIII	

ZIP	
	FRONT
	BACK

LEGEND: - MEASUREMENT UNAVAILABLE