


Custom measurement form:
PCO ORTHOSES CARD

PATIENT ID _____

PRODUCT CODE _____

NOTES _____




! WARNING!

- Chest and upper limb measurements should be taken with patient lying down.
- Hip and lower limb measurements are taken with patient standing.
- All measurements should be taken in centimeters!
- If the finger put between the child's body and the PCO product meets resistance, the size of the product has been selected correctly.

IMPORTANT!

Only a correctly filled in measurement card is a guarantee of proper compression of the product. We are not responsible for wrongly taken measurements!

HOW TO MEASURE?



- (I) Arms circumference at the level of the head of the humerus
- (II) Shoulder circumference measured through the armpit and acromion
- (III) Length measured from armpit to wrist
- (IV) Height measured from the armpit to the greater trochanter of the femur
- (V) Measurement from the manubrium of the sternum through the perineum to the cervical vertebra C7
- (VI) Height measured from waist to ankle
- (VII) Height measured from parineum to ankle
- (VIII) Patient is height

REQUIRED MEASUREMENTS

| PRODUCT | MEASUREMENT | LENGTH MEASUREMENT | CONTROL MEASUREMENT | ZIP |
|-------------------|--|------------------------|-------------------------|-----|
| HEAD | | | | |
| PCO-H-01 | 1 | | | |
| TRUNK | | | | |
| PCO-T-01 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18 | F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-02 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18 | F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-03 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18 | F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-04 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18 | F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-09 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18 | F, G, H, I, J | I, II, IV, VIII | Y |
| PCO-T-10 | 2, 3, 4, 5, 6, 7, 8, 14, 15, 16, 17, 18 | E, F, G, H, I, J | I, II, III, IV, VIII | Y |
| PCO-T-11 | 2, 3, 4, 5, 6, 7, 8, 9, 10, 14, 15, 16, 17, 18 | C, D, E, F, G, H, I, J | I, II, IV, VIII | Y |
| PCO-T-12 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19 | E, F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-13 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19 | F, G, H, I, J | I, II, III, IV, V, VIII | |
| PCO-T-14 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19, 20, 21 | F, G, H, I, J, K | I, II, IV, V, VII, VIII | |
| PCO-T-15 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19, 20, 21 | F, G, H, I, J, K | I, II, IV, V, VII, VIII | |
| PCO-T-16 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18 | F, G, H, I, J | I, II, IV, VIII | Y |
| PCO-T-17 | 2, 3, 4, 5, 6, 7, 8, 14, 15, 16, 17, 18 | E, F, G, H, I, J | I, II, III, IV, VIII | Y |
| PCO-T-18 | 2, 3, 4, 5, 6, 7, 8, 9, 10, 14, 15, 16, 17, 18 | C, D, E, F, G, H, I, J | I, II, III, IV, VIII | Y |
| PCO-T-19 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19 | F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-20 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19 | F, G, H, I, J | I, II, IV, V, VIII | |
| PCO-T-23 | 2, 3, 4, 5, 6, 14, 15, 16, 17, 18, 19 | F, G, H, I, J | I, II, IV, V, VIII | Y |
| PCO-T-24 | 2, 3, 4, 5, 6, 7, 8, 14, 15, 16, 17, 18, 19 | E, F, G, H, I, J | I, II, III, IV, V, VIII | Y |
| PCO-T-25 | 2, 3, 4, 5, 6, 7, 8, 9, 10, 14, 15, 16, 17, 18, 19 | C, D, E, F, G, H, I, J | I, II, III, IV, V, VIII | Y |
| UPPER LIMB | | | | |
| PCO-A-01 | 6, 7, 8, 9, 10, 11, 12, 13 | A, B, C, D, E | III | |
| PCO-A-02 | 8, 9, 10, 11, 12, 13 | A, B, C, D | III | |
| PCO-A-03 | 9, 10, 11, 12, 13 | A, B, C | III | |
| PCO-A-07 | 6, 7, 8, 9, 10 | C, D, E | III | |
| PCO-A-08 | 8, 9, 10 | C, D | III | |
| LOWER LIMB | | | | |
| PCO-L-01 | 16, 17, 18, 19, 20, 21, 22, 23, 24, 25 | I, J, K, L, M | VI, VII, VIII | |
| PCO-L-02 | 16, 17, 18, 19, 20, 21, 22, 23, 24 | I, J, K, L | VI, VII, VIII | |
| PCO-L-03 | 16, 17, 18, 19, 20, 21 | I, J, K | VI, VII, VIII | |
| PCO-L-08 | 19, 20, 21, 22, 23, 24, 25 | K, L, M | VII | |
| PCO-L-09 | 19, 20, 21, 22, 23, 24 | K, L | VII | |
| PCO-L-10 | 19, 20, 21 | K | VII | |
| PCO-L-11 | 23, 24, 25 | L, M | VII | |
| PCO-L-12 | 24, 25, 26, 27 | M, N, O, P | VII | |

| No. | MEASUREMENT | CIRCUMFERENCE | WIDTH | LENGTH |
|-----|---|---------------|-------|--------|
| | | | | |
| 1 | HEAD CIRCUMFERENCE | | | A |
| 2 | NECK CIRCUMFERENCE | | | B |
| 3 | NECK WIDTH | | | C |
| 4 | SHOULDERS WIDTH | | | D |
| 5 | SHOULDER-NECK WIDTH | | | E |
| 6 | ARM CIRCUMFERENCE | | | F |
| 7 | CIRCUMFERENCE BETWEEN 6 AND 8 | | | G |
| 8 | ELBOW CIRCUMFERENCE WITH A SLIGHTLY BENT HAND | | | H |
| 9 | CIRCUMFERENCE BETWEEN 8 AND 10 | | | I |
| 10 | WRIST CIRCUMFERENCE | | | J |
| 11 | THUMB CIRCUMFERENCE | | | K |
| 12 | HAND CIRCUMFERENCE | | | L |
| 13 | MCP JOINTS CIRCUMFERENCE | | | M |
| 14 | MEASUREMENT UNDER THE ARMPITS | | | N |
| 15 | MEASUREMENT BETWEEN 12 AND 14 | | | O |
| 16 | WAIST CIRCUMFERENCE | | | P |
| 17 | HIP CIRCUMFERENCE | | | |
| 18 | MEASUREMENT AT THE PERINEUM HEIGHT | | | |
| 19 | GROIN CIRCUMFERENCE | | | |
| 20 | CIRCUMFERENCE BETWEEN 19 AND 20 | | | |
| 21 | CIRCUMFERENCE 5 CM ABOVE THE KNEECAP | | | |
| 22 | KNEE CIRCUMFERENCE | | | |
| 23 | CIRCUMFERENCE 5 CM UNDER THE KNEECAP | | | |
| 24 | CALF CIRCUMFERENCE BETWEEN 22 AND 25 | | | |
| 25 | ANKLE CIRCUMFERENCE | | | |
| 26 | HEEL CIRCUMFERENCE | | | |
| 27 | METATARSUS CIRCUMFERENCE | | | |

| CONTROL MEASUREMENT | |
|---------------------|--|
| I | |
| II | |
| III | |
| IV | |
| V | |
| VI | |
| VII | |
| VIII | |

| ZIP | |
|-----|-------|
| | FRONT |
| | BACK |